


# BEACON Caresystem®

– a Physiological based Ventilator Assist and Monitoring System



The BEACON Caresystem® is a decision support system, implemented as an extra device on the side of the mechanical ventilator or as a stand-alone on a trolley, which provides advice as to appropriate ventilator settings. BEACON Caresystem® is based upon a number of linked mathematical models of physiology that enable simulation of changes in the patient on modifying ventilator settings.



## The BEACON Caresystem®

Learns about the individual patient's physiological status, which is reflected in model simulations and, consequently, the most appropriate ventilator settings for that individual, at that time.

# BEACON Caresystem®

A Physiological based Ventilator Assist and Monitoring System

## Physiological Based Ventilator Assist System

- The BEACON Caresystem® provides advice on the optimal  $\text{FiO}_2$ , P/Vt, Rf, I:E and PEEP setting regardless of patient clinical state
- The BEACON advice is based on Physiological Models, Clinical Preferences and Step2Target Simulations
- The BEACON Caresystem® is in contrast to rule-based systems, which use the same medical protocol or rules of thumb for all patients

## Intelligent Monitoring and Alarms

- The BEACON Caresystem's model parameters characterize the in-depth picture of the patient's state and physiology
- Monitoring is based on changes in physiology instead of a single point of measure, enabling improved patient safety
- Monitoring is in contrast to other systems, which measure and report a drop in  $\text{SpO}_2$ , but do not explain why the drop happened



Distributor label:




## A Physiological based Ventilator Assist and Monitoring System

MERMAID CARE

## Clinical challenges in Ventilation management



- Finding the optimal ventilator settings, based on the balance between physiological status and clinical preference
- Minimizing the risk of Ventilator-induced lung injury
- Weaning the patient from mechanical ventilation quickly and safely
- Balancing the goals of sufficient oxygenation and carbon dioxide removal



**MERMAID CARE**  
Mermaid Care A/S | Hedelund 1 | DK-9400 Nr. Sundby | Denmark | Tel. +45 7023 7015  
info@mermaidcare.com | www.mermaidcare.com | www.beaconcaresystem.com

MC.DOC.NO.00001805 #1.0



# Three different versions

The BEACON Caresystem® is available in three different versions with an easy upgrade path between them:

## BEACON D

Pulmonary diagnostic device for Intelligent Monitoring and Non-invasive V/Q measurement on ICU ventilated patients

## BEACON 3


Physiological Based Ventilation Assist and Intelligent Monitoring System, which provides advice for **FiO<sub>2</sub>, P/Vt, Rf** ventilator settings

## BEACON 5

Physiological Based Ventilation Assist and Intelligent Monitoring System, which provides advice for **FiO<sub>2</sub>, P/Vt, Rf, I:E and PEEP** ventilator settings

D

**DIAGNOSTIC OF SHUNT, LOW AND HIGH V/Q**  
**BEACON D:** Pulmonary diagnostic device for the Intelligent Monitoring and Non-invasive V/Q measurement on ICU ventilated patients




**Automatic Lung Parameter Estimation (ALPE)** provides an overview of pulmonary ventilation/perfusion (V/Q) distribution

- Pulmonary diagnostic, reporting SHUNT, Low V/Q and High V/Q distribution
- Historical review and trending for easy monitoring of changes in distribution
- Fast, non-invasive and guided measurement, using FiO<sub>2</sub> as tracer


V/Q distribution results can be displayed graphically or numerically making BEACON D an optimal Ventilator Assist tool

- Provides decision support for optimal ventilator settings
- Intuitive and easy-to-use interface
- Simple and minimal patient information required for measurement




**Physiological Based Ventilation Assist and Intelligent Monitoring System, which provides advice for three ventilator settings: FiO<sub>2</sub>, P/Vt and Rf**

- Intelligent Step2Target® algorithm for optimal FiO<sub>2</sub>, P/Vt and Rf setting advice
- Advice is based on the balance between patient's physiological status and clinical preferences
- Advice is provided for Volume/Pressure and Control/Support ventilation modes



**Intelligent Simulation, simulates the effect of ventilator settings**

- Simulates the effects of changes in FiO<sub>2</sub>, P/Vt and Rf
- Simulations describe the balance between
  - Over / under ventilation
  - stress / ventilator dependency
  - low oxygenation / oxygen toxicity
- Simulations presented on the PreferenceZone® hexagon, illustrating the balances
- Simulations provide a learning environment, to assess the effects of ventilation strategy



**Physiological Based Ventilation Assist and Intelligent Monitoring System, which provides advice for five ventilator setting: FiO<sub>2</sub>, P/Vt, Rf, I:E and PEEP**

- Intelligent Step2Target® algorithm for optimal FiO<sub>2</sub>, P/Vt, Rf, I:E and PEEP setting advice
- Advice is based on the balance between the patient's physiological status and clinical preferences
- Advice is provided for Volume/Pressure and Control/Support ventilation modes



**Intelligent Simulation, simulates the effect of ventilator settings**

- Simulates the effects of changes in FiO<sub>2</sub>, P/Vt and Rf, I:E and PEEP
- Simulations describe the balance between
  - Over / under ventilation
  - stress / ventilator dependency
  - low oxygenation / oxygen toxicity
- Simulations presented on the PreferenceZone® hexagon, illustrating the balances
- Simulations provide a learning environment, to assess the effects of ventilation strategy


3

**VENTILATOR ADVICE ON FIO<sub>2</sub>, P/Vt AND Rf**  
**BEACON 3:** Physiological Based Ventilation Assist and Intelligent Monitoring System connected to an ICU ventilator, which provides advice for three ventilator settings




**Intelligent Physiological Monitoring System**

- Monitoring of physiological changes (changes in V/Q distribution) and advice of when new ALPE measurement is needed
- Tracking and trending display of patient parameters
- Intuitive and easy-to-use interface with customizable layouts




**Intelligent Physiological Monitoring**

- Model parameters characterize the individual patient's physiological state
- Monitor changes based on physiology, and advises when a new ALPE measurement, Blood Gas Data or Cardiac Output data is needed, in order to tune BEACON's models to the patient's physiological condition
- Tracking and trending display of patient parameters
- Intuitive and easy-to-use interface with customizable layouts




**Intelligent Physiological Review**

- Tracking and trending of physiological changes, advice, patient parameters, ventilator settings etc.
- Intuitive and easy-to-use interface with customizable layouts
- Enables a complete overview of all past activities performed by the system



**Intelligent Physiological Monitoring System**

- Model parameters characterize the individual patient's physiological state
- Monitor changes based on physiology and advises when a new ALPE measurement, Blood Gas Data or Cardiac Output data is needed, in order to tune BEACON's models to the patient's physiological condition
- Tracking and trending display of patient parameters
- Intuitive and easy-to-use interface with customizable layouts



**Intelligent Physiological Review**

- Tracking and trending of physiological changes, advice, patient parameters, ventilator settings etc.
- Intuitive and easy-to-use interface with customizable layouts
- Enables a complete overview of all past activities performed by the system



Screen layout subject to changes without any further notice.  
Publications: See [www.beaconcaresystem.com/ beacon-d/publication](http://www.beaconcaresystem.com/ beacon-d/publication)

Screen layout subject to changes without any further notice.  
Publications: See [www.beaconcaresystem.com/ beacon-3/publication](http://www.beaconcaresystem.com/ beacon-3/publication)

Screen layout subject to changes without any further notice.  
Publications: See [www.beaconcaresystem.com/ beacon-5/publication](http://www.beaconcaresystem.com/ beacon-5/publication)